

Letter of Declaration

March 18, 2022

Dear customers,

This letter is served to notify you that COVID-19 Coronavirus Real Time PCR Kit, COVID-19 Coronavirus (ORF1ab/N/E) Real Time PCR Kit, COVID-19 Coronavirus and Influenza A/B Virus Real Time PCR Kit as confirmed can detect variants of SARS-CoV-2 [known as Alpha (B.1.1.7), Beta (B.1.351), Gamma (P.1), Delta (B.1.617.2, AY.1, AY.2, AY.3), Zeta (P.2), Eta (B.1.525), Iota (B.1.526), Kappa (B.1.617.1), Lambda (C.37), Fin-796H, N.9, B.1.616, Epsilon B.1.429, B.1.427, Theta (P.3), Mu (B.1.621), B.1.640, Omicron (B.1.1.529, BA.1, BA.1.1, BA.2, BA.3), IHU (B.1.640.2), BA.1xAY.4 recombinant (Delta-Omicron recombinant variant).

Targets of our related PCR kits are specific for ORF1ab gene, N gene, E gene. **NONE** of our kit detection targets are located at the reported mutation sites, which assures that the detection scale and performance will **NOT** be affected by these mutations.

We thank you for your continued support and understanding. As for questions; please kindly contact our sales representatives or send e-mail to info@biopertectus.com for more information.

JIANGSU BIOPERFECTUS TECHNOLOGIES CO., LTD

