

Letter of Declaration

May 16, 2022

Dear customers,

This letter is served to notify you that COVID-19 Coronavirus Real Time PCR Kit, COVID-19 Coronavirus (ORF1ab/N/E) Real Time PCR Kit, COVID-19 Coronavirus and Influenza A/B Virus Real Time PCR Kit as confirmed can detect variants of SARS-CoV-2 [known as Alpha (B.1.1.7), Beta (B.1.351), Gamma (P.1), Delta (B.1.617.2, AY.1, AY.2, AY.3), Zeta (P.2), Eta (B.1.525), Iota (B.1.526), Kappa (B.1.617.1), Lambda (C.37), Fin-796H, N.9, B.1.616, Epsilon B.1.429, B.1.427, Theta (P.3), Mu (B.1.621), B.1.640, IHU (B.1.640.2), Omicron (includes BA.1, BA.1.1, BA.2, BA.3, BA.4, BA.5 and descendent lineages), the recombinant lineage of BA.1 and BA.2, the recombinant lineage of Delta and Omicron].

Targets of our related PCR kits are specific for ORF1ab gene, N gene, E gene. **NONE** of our kit detection targets are located at the reported mutations, which assures that the detection scale and performance will **NOT** be affected by these mutations.

We thank you for your continued support and understanding. As for questions, please kindly contact our sales representatives or send e-mail to info@biopertectus.com for more information.

JIANGSU BIOPERPECTUS TECHNOLOGIES CO., LTD

